Dinbeat UNO

Protocolo de Telemedicina

INDICE

Protocolo Telemedicina

1. Consultas difíciles

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Beckman, B. (2013). Anesthesia and Pain Management for Small Animals. In *Veterinary Clinics of North America - Small Animal Practice* (Vol. 43, Issue 3, pp. 669–688). <https://doi.org/10.1016/j.cvsm.2013.02.006>

Blades Golubovic, S., & Rossmeisl Jr, J. H. (2017). Status epilepticus in dogs and cats, part 1: etiopathogenesis, epidemiology, and diagnosis. *Journal of Veterinary Emergency and Critical Care*, 27(3), 278–287.

Blades Golubovic, S., & Rossmeisl, J. H. (2017). Status epilepticus in dogs and cats, part 2: treatment, monitoring, and prognosis. *Journal of Veterinary Emergency and Critical Care*, 27(3), 288–300. <https://doi.org/10.1111/vec.12604>

Brambilla, P. G., Polli, M., Pradelli, D., Papa, M., Rizzi, R., Bagardi, M., & Bussadori, C. (2020). Epidemiological study of congenital heart diseases in dogs: Prevalence, popularity, and volatility throughout twenty years of clinical practice. *PLoS ONE*, 15(7 July). <https://doi.org/10.1371/journal.pone.0230160>

Bruno, E. A., Guthrie, J. W., Ellwood, S. A., Mellanby, R. J., & Clements, D. N. (2015). Global positioning System derived performance measures are responsive indicators of physical activity, disease, and the success of clinical treatments in domestic dogs. *PLoS ONE*, 10(2). <https://doi.org/10.1371/journal.pone.0117094>

Camps, T., Amat, M., & Manteca, X. (2019). A review of medical conditions and behavioral problems in dogs and cats. In *Animals* (Vol. 9, Issue 12). MDPI AG. <https://doi.org/10.3390/ani9121133>

Chalifoux, N. v., Drobatz, K. J., & Reineke, E. L. (2021). Predictors of inflammatory lower airway disease in cats presented to the emergency room in respiratory distress: a case-control study. *Journal of Feline Medicine and Surgery*, 23(12), 1098–1108. <https://doi.org/10.1177/1098612X21996145>

Dijkstra, E., Teske, E., & Szatmári, V. (2018). Respiratory rate of clinically healthy cats measured in veterinary consultation rooms. *Veterinary Journal*, 234, 96–101. <https://doi.org/10.1016/j.tvjl.2018.02.014>

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Downing s, R. (2011). Pain management for veterinary palliative care and hospice patients. In *Veterinary Clinics of North America - Small Animal Practice* (Vol. 41, Issue 3, pp. 531–550). <https://doi.org/10.1016/j.cvsm.2011.03.010>

Dyson, D. H. (2008). Perioperative Pain Management in Veterinary Patients. In *Veterinary Clinics of North America - Small Animal Practice* (Vol. 38, Issue 6, pp. 1309–1327). <https://doi.org/10.1016/j.cvsm.2008.06.006>

Fan, T. M. (2014). Pain management in veterinary patients with cancer. In *The Veterinary clinics of North America. Small animal practice* (Vol. 44, Issue 5, pp. 989–1001). <https://doi.org/10.1016/j.cvsm.2014.05.005>

Fuentes, L., Johnson, L. R., & Dennis, S. (2010). *BSAVA Manual of canine and feline cardiorespiratory medicine*. British Small Animal Veterinary Association.

Harper, T. A. M. (2017). Conservative Management of Hip Dysplasia. In *Veterinary Clinics of North America - Small Animal Practice* (Vol. 47, Issue 4, pp. 807–821). W.B. Saunders. <https://doi.org/10.1016/j.cvsm.2017.02.007>

Höglund, K., Hanås, S., Carnabuci, C., Ljungvall, I., Tidholm, A., & Häggström, J. (2012). Blood Pressure, Heart Rate, and Urinary Catecholamines in Healthy Dogs Subjected to Different Clinical Settings. *Journal of Veterinary Internal Medicine*, 26(6), 1300–1308. <https://doi.org/10.1111/j.1939-1676.2012.00999.x>

Ljungvall, I., Rishniw, M., Porciello, F., Häggström, J., & Ohad, D. (2014). Sleeping and resting respiratory rates in healthy adult cats and cats with subclinical heart disease. *Journal of Feline Medicine and Surgery*, 16(4), 281–290. <https://doi.org/10.1177/1098612X13508940>

Loewen, J. M., & Bach, J. F. (2022). Respiratory distress in small animals: Pathophysiology and clinical approach. *Journal of Veterinary Emergency and Critical Care*, 32, 3–15. <https://doi.org/10.1111/vec.13121>

Looney, A. (2010). Oncology Pain in Veterinary Patients. In *Topics in Companion Animal Medicine* (Vol. 25, Issue 1, pp. 32–44). <https://doi.org/10.1053/j.tcam.2009.10.008>

Mandese, W. W., Griffin, F. C., Reynolds, P. S., Blew, A. C., Deriberprey, A. S., & Estrada, A. H. (2021). Stress in client-owned dogs related to clinical exam location: a randomised crossover trial. *Journal of Small Animal Practice*, 62(2), 82–88. <https://doi.org/10.1111/jsap.13248>

Moore, S. A. (2016). Managing neuropathic pain in dogs. In *Frontiers in Veterinary Science* (Vol. 3, Issue FEB). Frontiers Media S.A. <https://doi.org/10.3389/fvets.2016.00012>

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Petrie, J. P. (2005). Practical application of holter monitoring in dogs and cats. *Clinical Techniques in Small Animal Practice*, 20(3 SPEC. ISS.), 173–181. <https://doi.org/10.1053/j.ctsap.2005.05.006>

Porciello, F., Rishniw, M., Ljungvall, I., Ferasin, L., Haggstrom, J., & Ohad, D. G. (2016). Sleeping and resting respiratory rates in dogs and cats with medically-controlled left-sided congestive heart failure. *Veterinary Journal*, 207, 164–168. <https://doi.org/10.1016/j.tvjl.2015.08.017>

Rae, L., MacNab, N., Bidner, S., Davidson, C., & McDonagh, P. (2021). Attitudes and practices of veterinarians in Australia to acute pain management in cats. *Journal of Feline Medicine and Surgery*. <https://doi.org/10.1177/1098612X211043086>

Reinero, C., Visser, L. C., Kellihan, H. B., Masseau, I., Rozanski, E., Clercx, C., ... & Scansen, B. A. (2020). ACVIM consensus statement guidelines for the diagnosis, classification, treatment, and monitoring of pulmonary hypertension in dogs. *Journal of veterinary internal medicine*, 34(2), 549-573.

Rishniw, M., Ljungvall, I., Porciello, F., Häggström, J., & Ohad, D. G. (2012). Sleeping respiratory rates in apparently healthy adult dogs. *Research in Veterinary Science*, 93(2), 965–969. <https://doi.org/10.1016/j.rvsc.2011.12.014>

Sigrist, N. E., Adamik, K. N., Doherr, M. G., & Spreng, D. E. (2011). Evaluation of respiratory parameters at presentation as clinical indicators of the respiratory localization in dogs and cats with respiratory distress. *Journal of Veterinary Emergency and Critical Care*, 21(1), 13–23. <https://doi.org/10.1111/j.1476-4431.2010.00589.x>

Srithunyarat, T., Höglund, O. v., Hagman, R., Olsson, U., Stridsberg, M., Lagerstedt, A. S., & Pettersson, A. (2016). Catestatin, vasostatin, cortisol, temperature, heart rate, respiratory rate, scores of the short form of the Glasgow composite measure pain scale and visual analog scale for stress and pain behavior in dogs before and after ovariohysterectomy. *BMC Research Notes*, 9(1). <https://doi.org/10.1186/s13104-016-2193-1>

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Steagall, P. v., Robertson, S., Simon, B., Warne, L. N., Shilo-Benjamini, Y., & Taylor, S. (2022). 2022 ISFM Consensus Guidelines on the Management of Acute Pain in Cats. *Journal of Feline Medicine and Surgery*, 24(1), 4–30. <https://doi.org/10.1177/1098612X211066268>

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